

Welch®
Laboratory Vacuum Catalog
2005-2006

Typical Set-ups for Lab Applications
Catalog Page 4

ORDER ONLINE AT: www.cfpwarehouse.com or contact:

COMBINED FLUID PRODUCTS COMPANY
805 Oakwood Road, Lake Zurich, Illinois 60047
sales@combinedfluidproducts.com
www.combinedfluidproducts.com
Tel (847) 540-0054 Fax (847) 540-0513
(800) 521-2083

Typical Set-ups For Lab Applications

Proper vacuum system set-up in lab applications reduces pump maintenance which lowers cost and increases user productivity. Proper vacuum system set-ups consist of the lab appliance, tubing, vacuum pump, and normally a cold or liquid trap. Failure to use or maintain the appropriate trap will result in a shortened service life of the pump (See page 32 for trap options).

